Uncertainties for coincidence counting

Siobhan Tobin

November 10, 2017

The following convention is used:

- $\gamma = \text{gamma}$, $\beta = \text{ABC}$ or LD, C = ABCX or LDX, depending on whether doubles or triples data is being used;
- C refers to raw counts, e.g. C_{β} is the number of beta counts;
- R refers to uncorrected rate, e.g. R_{γ} is the gamma rate;
- B refers to background rate, e.g. B_C is the background coincidence rate;
- N refers to corrected rate, e.g. N_C is the decay- and background-corrected coincidence count rate.

A is the dilution factor, D is the decay factor, and m is the mass of the source. Negligible uncertainty in A, D, and m is assumed throughout.

β count rate per unit mass

The β count rate per unit mass (here referred to as F) is

$$F = N_{\beta} \cdot \frac{A}{m}$$

with uncertainty

$$u_F = \sqrt{\frac{1}{C_{\beta}} + \frac{1 - N_{\beta}/N_{\gamma}}{C_C}} \cdot N_{\beta} \cdot \frac{A}{m}.$$

In [?], the following formula is given:

$$\frac{u_F}{F} = \left(\frac{1 - \epsilon_\beta}{N_C}\right)^{1/2}$$

where ϵ_{β} is the β detection efficiency.

$\beta \gamma/C$ count rate per unit mass

The $\beta \gamma / C$ count rate per unit mass (here referred to as H) is

$$H = \frac{N_{\beta} N_{\gamma}}{N_C} \cdot \frac{A}{m}.$$

The uncertainty is

$$u_H = \sqrt{\frac{1}{C_\beta} + \frac{1 - N_\beta/N_\gamma}{C_C}} \cdot \frac{N_\beta N_\gamma}{N_C} \cdot \frac{A}{m}.$$

In [?], the following formula is given:

$$\frac{u_H}{H} = \left(\frac{1}{n_C}(2\epsilon_\beta\epsilon_\gamma - \epsilon_\beta - \epsilon_\gamma + 1)\right)^{1/2}$$

where ϵ_{γ} is the γ detection efficiency.

$\gamma/C-1$ "inverse efficiency" parameter

Here $\gamma/C - 1$ is defined as

$$J = \frac{N_{\gamma}}{N_C} - 1$$

The uncertainty is complicated by the uncertainty in the background counts.

$$u_J^2 = \left(\frac{\partial J}{\partial N_\gamma}\right)^2 u_{N_\gamma}^2 + \left(\frac{\partial J}{\partial N_C}\right)^2 u_{N_C}^2 \tag{1}$$

$$= \frac{1}{N_C^2} u_{N_\gamma}^2 + \frac{N_\gamma^2}{N_C^4} u_{N_C}^2 \tag{2}$$

The γ count rate is given by

$$N_{\gamma} = \left(\frac{C_{\gamma}}{\tau} - B_{\gamma}\right) \cdot D$$

The uncertainty in the background γ count rate is simply

$$u_{B_{\gamma}} = \sigma_{B_{\gamma}}$$
.

Similarly,

$$u_{B_C} = \sigma_{B_C}$$
.

$1-C/\gamma$ "inefficiency" parameter

Here $1 - C/\gamma$ is defined as

$$K = 1 - \frac{N_C}{N_{\gamma}}.$$

The uncertainty is

$$u_K^2 = \left(\frac{\partial K}{\partial N_\gamma}\right)^2 u_{N_\gamma}^2 + \left(\frac{\partial K}{\partial N_C}\right)^2 u_{N_C}^2 \tag{3}$$

$$=\frac{N_C^2}{N_\gamma^4}u_{N_\gamma}^2 + \frac{1}{N_\gamma^2}u_{N_C}^2 \tag{4}$$

References

[1] P. J. Campion and J. G. V. Taylor, Statistical Errors in Disintegration Rate Measurements by the Coincidence Technique, Int. J. Appl. Radiat. Isot. 10 131-133 (1961)